



## DK50 2V/50

For two dental units  
 Pour deux unités dentaires  
 Para dos unidades dentales



## DK50 2x2V/110

For four dental units  
 Pour quatre unités dentaires  
 Para cuatro unidades dentales



Accessories/ Accessoires/ Accesorios: **BOX MD KJF FS REG AOK**

	DK50 2V/50	DK50 2V/50/M	DK50 2V/50S	DK50 2V/50S/M
*1	2	2	2	2
DRYER	-	✓	-	✓
BOX	-	-	✓	✓
(5-8bar)	135 l.min <sup>-1</sup> / 6bar	104 l.min <sup>-1</sup> / 6bar	135 l.min <sup>-1</sup> / 6bar	104 l.min <sup>-1</sup> / 6bar
(8-10bar)	110 l.min <sup>-1</sup> / 8bar	90 l.min <sup>-1</sup> / 8bar	110 l.min <sup>-1</sup> / 8bar	90 l.min <sup>-1</sup> / 8bar
S1-S3	S1 - 100%	S1 - 100%	S1 - 100%	S1 - 100%
	≤ 67 dB(A)	≤ 69 dB(A)	≤ 51,5 dB(A)	≤ 55 dB(A)

Accessories/ Accessoires/ Accesorios: **BOX MD KJF FS REG AOK**

	DK50 2x2V/110	DK50 2x2V/110/M	DK50 2x2V/110S	DK50 2x2V/110S/M
*1	4	3 - 4	4	3 - 4
DRYER	-	✓	-	✓
BOX	-	-	✓	✓
(5-7bar)	280 l.min <sup>-1</sup> / 5bar	225 l.min <sup>-1</sup> / 5bar	280 l.min <sup>-1</sup> / 5bar	225 l.min <sup>-1</sup> / 5bar
(6-8bar)	270 l.min <sup>-1</sup> / 6bar	215 l.min <sup>-1</sup> / 6bar	270 l.min <sup>-1</sup> / 6bar	215 l.min <sup>-1</sup> / 6bar
(8-10bar)	210 l.min <sup>-1</sup> / 8bar	165 l.min <sup>-1</sup> / 8bar	210 l.min <sup>-1</sup> / 8bar	165 l.min <sup>-1</sup> / 8bar
S1-S3	S1 - 100%	S1 - 100%	S1 - 100%	S1 - 100%
	≤ 73 dB(A)	≤ 73 dB(A)	≤ 55 dB(A)	≤ 56 dB(A)