

# Triton Stool

*with backrest*



You can easily adjust the height of the backrest.



The castors are equipped with a special rubber lining to move smoothly without straining the floor surfaces.

The seat has a metallic, smooth surface, a base for easy cleaning, and a high-strength gas spring for adjusting the desired height.

The upholstery is molded, seamless, in a large stock of spare parts. This offers the advantage of easy and extremely economical replacement in case of accidental wear.

All-metal construction compact and robust, specially designed for professional use.

Soft foam pillow rich in material.



Available in 2 sizes : Short and Medium

| Size   | Minimum height* | Maximum height* |
|--------|-----------------|-----------------|
| SHORT  | 47 cm           | 61 cm           |
| MEDIUM | 56 cm           | 76 cm           |

\*\* The floor surface distance from the seat surface in centimeters

## COLORS AVAILABLE



\* Metallic Shades