





# OWOODY-RXPORTABLE

LESS WEIGHT FOR BETTER IMAGES
PORTABLE DENTAL X-RAY SYSTEM

# OWANDY-RXPORTABLE





#### **USER FRIENDLY**

Parameters are displayed in a simple way **to facilitate use**: battery level, type of examination and patient, X-ray emission, kV, mA... You can also configure **Owandy-RX Portable** according to your needs thanks to the "MEMO" function.





#### **LIGHTWEIGHT**

Weighing 1.8kg, **Owandy-RX Portable** is one of the lightest X-ray devices on the market. It is therefore easy to handle and to position.



#### LONG BATTERY LIFE

With one charge, you can shoot up to 300 images!



#### **HIGH DEFINITION IMAGES**

Thanks to its **0.3mm focal spot**, you will obtain high quality images, both with sensor and plates.









### **Technical** features

X-ray tube	CEI
kV/mA	70kV / 3mA
Focal spot	0.3 mm
Total filtration	1.5mm Al
Battery charger input voltage	100-240V / 50~60Hz
Battery voltage	24.2V
Weight	1.8 kg
Options	Backscatter protection Rectangular collimator (30x40mm, 20x30 mm)

## DIGITAL WORKFLOW OWANDY RADIOLOGY

A COMPREHENSIVE RANGE TO MEET ALL YOUR REQUIREMENTS



